

SWIFT OPTICAL INSTRUMENTS, INC.

Microscopes • Digital Imaging Products



"Quality Microscopes for a Lifetime"

M15 Series

Infinity Corrected Microscope

Performance

- ◆ Swift high-resolution infinity corrected optics are parfocaled and parcentered
- ◆ Slow-closing, slide holder clip on mechanical stage minimizes slide breakage
- ◆ Plastic iris diaphragm shield protects iris from damage
- ◆ 20mm eyepieces for increased field of view

Precision

- ◆ Siedentopf Head provides larger interpupillary range and is ergonomically positioned for operator comfort
- ◆ Focus-tension control knob offers "feather-touch" sensitivity for extremely precise focusing
- ◆ High-grade synthetic lubricants provide a smooth, precise operation

Ease of Use

- ◆ Heads are inclined at 45° for ergonomic viewing
- ◆ Inward-facing (reversed) nosepiece allows user clear access to stage for slide manipulation
- ◆ 40XRD & 100XRD objectives are spring-loaded and sealed to protect optics from contact or oil damage
- ◆ LED illumination with variable intensity provides a energy-efficient, cool, white light
- ◆ Abbe 1.25 NA brightfield condenser is optically coated for maximum light transmission and enhanced contrast

Quality

- ◆ Swift quality construction assures "trouble free" operation and minimal down time
- ◆ All-metal internal parts for durability
- ◆ Swift Optical Limited Lifetime Warranty covers defective parts and workmanship for the life of the instrument



M15 Series

MODEL	HEAD	EYEPIECE	OBJECTIVE	CONDENSER	ILLUMINATOR
M15	Binocular	W10XD, 20mm	4XD, 10XD, 40XRD, 100XRD Infinity corrected	Abbe 1.25 DIN Rack & Pinion	Variable LED

Modern Finish: Unless otherwise indicated, all instruments are finished in acid and reagent resistant, Swift epoxy-ester resin.

The Swift Optical Instruments, Inc. Limited Lifetime Warranty assures that the microscope is guaranteed against defects in material and workmanship for the life of the product. Electrical components are covered for 5 year after purchase. Video and camera components are covered for 1 year after purchase. Normal wear, routine maintenance, light bulbs or damage resulting from repair by unauthorized parties, accident, alteration, shipping, misuse or abuse is not covered. Warranty service is provided by Swift Optical Instruments, Inc.'s authorized technicians. Determination of warranty is at the technician's discretion. Defective products covered by the warranty will be repaired free of charge when they are returned, postpaid, to Swift Optical headquarters. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. Limited to USA.

Service: Swift Optical maintains a highly-skilled and specially-trained team of technicians who are available to make any needed repairs or adjustments.

Design Changes: Swift Optical reserves the right to make changes in instrument design, in accordance with scientific and mechanical progress, without notice or obligation.

NOTES:



www.Swift-MicroscopeWorld.com
800-942-0528 Toll Free
760-438-0528 International
info@swift-microscopeworld.com

SWIFT OPTICAL INSTRUMENTS, INC.

Customer Service: 877.967.9438

www.swiftoptical.com

