

# SWIFT OPTICAL INSTRUMENTS, INC.

Microscopes • Digital Imaging Products

"Quality Microscopes for a Lifetime"



## M3 Series

### Micro/ Macro/ Field Microscopes

#### Versatility

- ◆ Three microscopes in one. Can be used as a microscope, macroscope or in the field.
- ◆ Lends itself to a wide variety of field activities from gross dissection to cell identification
- ◆ Re-chargable illumination allows up to 40 hours of use on a single 8 hour charge

#### Unique Design

- ◆ Sturdy detachable tri-pod legs fold for easy storage
- ◆ Multiple stage positions allow for micro and macro viewing
- ◆ Includes stage plate for viewing micro specimens and black/white contrast plate and clear collection container for viewing macro specimens
- ◆ Ergonomic carry handle, promotes proper handling

#### LED Illumination

- ◆ Top and bottom LED illumination provides energy efficient, white light with virtually no heat
- ◆ Extra bright LED lasts up to 50,000 hours between replacements for minimum maintenance and downtime

#### Quality

- ◆ Quality manufacturing and all-metal internal parts ensure durability and minimal down time
- ◆ Swift Optical Limited Lifetime Warranty covers defective parts and workmanship for the life of the instrument



M3-B



M3-M

NEW MODELS!



[www.Swift-microscopeworld.com](http://www.Swift-microscopeworld.com)  
800-942-0528 Toll Free  
760-438-0528 International  
[info@swift-microscopeworld.com](mailto:info@swift-microscopeworld.com)

SWIFT OPTICAL INSTRUMENTS, INC.

Customer Service: 877.967.9438  
[www.swiftoptical.com](http://www.swiftoptical.com)

# Swift M3 Series

MODEL	HEAD	EYEPIECE	OBJECTIVE	STAGE	ILLUMINATOR
M3-B	Binocular	W10XD, 18mm	4XD, 10XD, 40XRD Micro, 2X Macro	Stage Plate, Specimen Cup, Contrast Plate	Cordless, Top and Bottom LED
M3-M	Monocular	W10XD, 18mm	4XD, 10XD, 40XRD Micro, 2X Macro	Stage Plate, Specimen Cup, Contrast Plate	Cordless, Top and Bottom LED

Unique adjustable stage position allows micro and macro viewing



### Micro mode

Align stage with micro stage position dot and use the stage plate for viewing slides



Transitioning from micro/macro mode is as easy as loosening the stage thumbscrew, re-positioning the stage, re-tightening the thumbscrew, and then swapping stage inserts.



### Macro mode

Align stage with a macro stage position dot and use the contrast plate or collection cup to view macro specimens

**Modern Finish:** Unless otherwise indicated, all instruments are finished in acid and reagent resistant, Swift epoxy-ester resin.

**Limited Lifetime Warranty:** The Swift Optical Instruments, Inc. Limited Lifetime Warranty assures that the microscope is guaranteed against defects in material and workmanship for the life of the product. Electrical components are covered for three years; camera and video components are covered for one year after purchase. Normal wear, routine maintenance, light bulbs or damage resulting from repair by unauthorized parties, accident, alteration, shipping, misuse or abuse is not covered. Warranty service is provided by Swift Optical Instruments, Inc. authorized technicians. Determination of warranty is at the technician's discretion.

Other than set forth above, Swift hereby disclaims all warranties, expressed or implied, of fitness for a particular purpose.

**Service:** Swift Optical maintains a highly-skilled and specially-trained team of technicians who are available to make any needed repairs or adjustments.

**Design Changes:** Swift Optical reserves the right to make changes in instrument design, in accordance with scientific and mechanical progress, without notice or obligation.



[www.Swift-Microscopeworld.com](http://www.Swift-Microscopeworld.com)  
800-942-0528 Toll Free  
760-438-0528 International  
[info@swift-microscopeworld.com](mailto:info@swift-microscopeworld.com)



**SWIFT OPTICAL INSTRUMENTS, INC.**

Customer Service: 877.967.9438  
[www.swiftoptical.com](http://www.swiftoptical.com)